

EXACT EXCHANGER I N C O R P O R A T E D

Sept. 28 2020

Tescorp
PO Box 471453
Tulsa, OK 74147-1453

Dana,

For your record, all of the following SN's were tested successfully with dry Nitrogen at 1.1 x the design pressure of 25 PSIG and were held under pressure with no leaks for a minimum of 30 minutes.

Tescorp PO Number: F-0820-0351
Job Number: 320015
Serial Numbers: 2009656

Should you have any questions please let us know.

Regards,
David Self